



## Preparing the Medium for Use and Storage

- The medium and supplements are delivered cooled (@ 2 to 8°C).
- Store the medium @ 4 to 8°C in the dark.
- The supplements should be frozen @ -20°C after arrival.
- It is recommended to aliquot the supplements before and keep them frozen @ -20°C until use.
- If stored properly, the products are stable until the expiry date.
- To give the cultivation media thaw the supplements at 15°C to 25°C and apply the solution to the basal medium (it is recommended to prepare small batches of complete medium).
- The complete medium is stable for approximately one month when stored at 4°C.